

MIqPCR: *Minimum Information about a Quantitative Polymerase Chain Reaction experiment*

Version 1.0, 4th April, 2008.

RDML-Consortium

This module identifies the minimum information required to report the use of quantitative Polymerase Chain Reaction (qPCR) assay in a manner sufficient to support the unambiguous interpretation of the results presented and the potential corroboration of the conclusions drawn.

Introduction

Quantitative polymerase chain reaction (qPCR) assays measure the copies of a specific DNA target in a sample as that sample is repeatedly passed through the polymerase chain reaction. Special qPCR machines are required to quantify the amplification products at each step of the cycle. MIqPCR describes the minimal information needed to allow effective interpretation and reanalysis of the data generated by such assays. For further background on this and related projects appropriate for qPCR, please visit the website† (<http://www.rdml.org>).

The following section, detailing the reporting requirements for the use of qPCR, is subdivided as follows:

1. Administrative information
2. Sample annotation
3. Target annotation
4. Thermal Cycling Conditions
5. Run data
6. Software requirements

Reporting requirement for Quantitative PCR Assays

1. Administrative information

- a) Experiment description
 - Experiment description
 - Responsible person and contact details

2. Sample annotation

- a) Sample description
 - Sample ID
 - Sample description
 - cDNA synthesis method and DNase treatment (cDNA samples only)
 - Template quantity (standard and optical calibrator samples only)
- b) Sample role in qPCR assay
 - Sample type
 - Inter run calibrator (true or false)
 - Calibrator sample (true or false)

3. Target annotation

- a) Target description
 - Target ID

- Sequence of primers OR commercial assay description
- b) Target role in qPCR assay
 - Target type

4. Thermal Cycling Conditions Information

- a) PCR program
 - Complete description of the cycling conditions

5. Run data

- a) Instrument information
 - Plate format
 - Instrument description
 - Software description and version
- b) Information required for each well
 - Well ID
 - Sample ID
 - Target ID
 - Amplification curve fluorescence values for each data point
 - Melting curve fluorescence values for each data point
 - Quantification Cycle

6. Software requirements

- a) RDML-Support
 - Software solutions, including databases, must support the import and export of RDML files.
 - qPCR machines must allow the export of raw data for the amplification as well as for melting curves.

Summary

These minimum reporting requirements for quantitative PCR assays specify that a significant degree of detail be captured, for sample and target description, PCR amplification and data capture and analysis, and that in addition the raw and processed data are also made available. However, it is clear that providing the information required by this document will enable the effective interpretation and assessment of qPCR data and protocols and potentially, support experimental corroboration. Much of the information required herein may already be stored in an electronic format, or exportable from the instrument; we anticipate further automation of this process.

These guidelines will evolve. To contribute, or to track the process to remain 'MIqPCR-compliant', browse to the website at <http://www.rdml.org/>

Corresponding author

Jan Hellemans, PhD
Center for Medical Genetics Ghent (CMGG)
Ghent University Hospital
Medical Research Building (MRB), 1st floor, room 110.050
De Pintelaan 185, B-9000 Ghent, Belgium
Phone : +32 9 33 25535
Fax : +32 9 33 26549
Jan.Hellemans@UGent.be

Glossary

- Experiment description: A sentence or short phrase describing the context of the experiment, or a stable pointer to contextualizing information, such as a database record or a publication.
- Responsible person and contact details: Contact details of the experiment supervisor, this includes the name and email address, and name and postal address of the institute
- Sample ID: Unique identifier for a given sample, to be used in the run data as reference for the sample.
- Sample description: A sentence or short phrase describing the sample.
- cDNA synthesis method: Enzyme and priming method used for cDNA synthesis. Possible priming methods are: *oligo-dt*, *random*, *target specific* or *any combination of the above*.
- DNase treatment: Boolean value to indicate whether nucleic acids have been treated or not with DNase. Possible values are *true* and *false*.
- Template quantity: Nucleic acid quantity of a standard sample
- Sample type: Describes the type of sample. Possible values are : *unknown*, *no template control*, *no amplification control*, *standard* or *optical calibrator*.
- Inter-run calibrator: Boolean value to indicate whether or not a sample can be used to detect and correct inter-run differences. Possible values are *true* and *false*.
- Calibrator sample : Indicates if a sample is used as a reference to rescale all relative measurements to (e.g. untreated control). Possible values are *true* and *false*.
- Target ID: Unique identifier for a given target, to be used in the run data as a reference to the detailed target information.
- Sequence of primers: Nucleotide sequence of the PCR primers. The amplicon sequence and optional 3'- and 5'-tags can also be mentioned.
- Commercial assay description: Short description of the commercial assay used for PCR amplification, includes at least company name and order number.
- Target type: Describes the use of a target. Possible values are: *reference* and *target of interest*.
- PCR program: Describes the thermal profile used for the PCR amplification and detection of accumulating amplicons. The program consist of one or multiple steps, each with it own characteristic information (e.g. duration, temperature, type of measurement, a loop with number of repeats, etc.)
- PCR format: Description of the run format used in the experiment (for example 96-well plate A1-H12, or 32-well rotor 1-32).
- Instrument description: Name of instrument and vendor
- Software description and version: Contains the name and version of the software used for the analysis. The methods used for determining the background and calculating the quantification cycle can be included as well.

Well ID: Identifier for a given well. Format is free, but should preferably reflect the run format.

Amplification curve fluorescence values: Contains the amplification fluorescence data for a given amplicon. The data is organized per sample and per amplicon with a corresponding cycle number for each data point.

Melting curve fluorescence values: Contains the fluorescence data of the melting curve for a given amplicon. The data is organized per sample and per amplicon with a corresponding temperature for each data point.

Quantification Cycle (Cq): Universal name for the fractional PCR cycle at which the target is quantified in a given sample. Depending on the real-time instrument, either crossing point (Cp), cycle threshold (Ct) or a take-off point (Top) are used to refer to the same value (Cq).

Intellectual Property Statement

The RDML-Consortium takes no position regarding the validity or scope of any intellectual property or other rights that might be claimed to pertain to the implementation or use of the technology described in this document or the extent to which any license under such rights might or might not be available; neither does it represent that it has made any effort to identify any such rights. Copies of claims of rights made available for publication and any assurances of licenses to be made available, or the result of an attempt made to obtain a general license or permission for the use of such proprietary rights by implementers or users of this specification can be obtained from the RDML-Consortium.

The RDML-Consortium invites any interested party to bring to its attention any copyrights, patents or patent applications, or other proprietary rights which may cover technology that may be required to practice this recommendation. Please address the information to the RDML-Consortium (see contacts information at RDML-Consortium website).

Copyright Notice

Copyright (C) RDML-Consortium (2008). All Rights Reserved.

This document and translations of it may be copied and furnished to others, and derivative works that comment on or otherwise explain it or assist in its implementation may be prepared, copied, published and distributed, in whole or in part, without restriction of any kind, provided that the above copyright notice and this paragraph are included on all such copies and derivative works. However, this document itself may not be modified in any way, such as by removing the copyright notice or references to the RDML-Consortium or other organizations, except as needed for the purpose of developing qPCR Recommendations in which case the procedures for copyrights defined in the RDML-Consortium Document process must be followed, or as required to translate it into languages other than English.

The limited permissions granted above are perpetual and will not be revoked by the RDML-Consortium or its successors or assigns.

This document and the information contained herein is provided on an "AS IS" basis and THE RDML CONSORTIUM DISCLAIMS ALL WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO ANY WARRANTY THAT THE USE OF THE INFORMATION HEREIN WILL NOT INFRINGE ANY RIGHTS OR ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE."

